

CapMed+

1995 OEC 9600 ESP



1995 OEC C-Arm ESP with 9"ii

Unit is staged and available for viewing at our facility

- Serial Number 62-1495
- Good working condition
- New overlays
- New power supply and cord
- YOM 1995

System details

OEC 9600 ESP 9" image intensifier.

50 Hz camera providing 12 fps

Software version and revision

- Control panel V1.24 Rev 8
- W/S Station 876780 Rev 18

Cine Disk: Alta 15

Line frequency: 50hz

- Available for pick up at our facility in Tulsa, OK. Shipping and installation is available.
- Sold as-is where-is. No warranty

Actual Photos





Screen Shots

